

	<b>Package contents</b>	<b>Catalog Number</b> 18091050, 18091150 18091200, 18091300	<b>Size</b> 50 reactions 200 reactions	Kit Contents
	<b>Storage conditions</b>	Store at -20°C (non-frost-free)		
	<b>Required materials</b>	<ul style="list-style-type: none"> <li>▪ Template: RNA</li> <li>▪ <i>(Optional)</i> 2 µM gene-specific primers</li> </ul>		
	<b>Timing</b>	<ul style="list-style-type: none"> <li>▪ Preparation time: 10 minutes</li> <li>▪ Run time: 20 minutes</li> </ul>		
	<b>Selection guides</b>	Go online to learn more. <a href="#">PCR Enzymes and Master Mixes</a> <a href="#">RT Enzymes and Kits</a> <a href="#">Real-Time PCR Instruments</a> <a href="#">Real-Time PCR Master Mixes</a> <a href="#">PCR Thermal Cyclers</a>		
	<b>Product description</b>	For first strand cDNA synthesis using total RNA or poly(A)+-selected RNA primed with oligo(dT), random primers, or a gene-specific primer.		
	<b>Important guidelines</b>	Pre-warm the 5× SSIV Buffer to room temperature before use. Vortex and briefly centrifuge the buffer prior to preparing the reverse transcription reaction mix.		
	<b>Online resources</b>	Visit our <a href="#">product page</a> for additional information and protocols. For support, visit <a href="http://thermofisher.com/support">thermofisher.com/support</a> .		

## Protocol outline

- A. *(Optional)* Remove genomic DNA
- B. Anneal primers to RNA
- C. Assemble reaction mix
- D. Add reaction mix to annealed RNA

## RT reaction setup

Use the measurements below to prepare your RT reaction, or enter your own parameters in the column provided.

Component	20-µL rxn	Custom	Final Conc.
DEPC-treated water	to 20 µL	to µL	N/A
5× SSIV Buffer	4.0 µL	µL	1×
10 mM dNTP mix (10 mM each)	1.0 µL	µL	0.5 mM each
100 mM DTT	1.0 µL	µL	5 mM
Ribonuclease Inhibitor (40 U/µL)	1.0 µL	µL	2.0 U/µL
50 µM Oligo d(T) <sub>20</sub> primer, or 50 ng/µL random hexamers, or 2 µM gene-specific primer	1.0 µL 1.0 µL 1.0 µL	µL	2.5 µM 2.5 ng/µL 0.1 µM
Template RNA*	varies	µL	< 5 µg total RNA or < 500 ng mRNA

\* 10 pg–5 µg total RNA or 10 pg–500 ng mRNA

## RT protocol

- Go to page 3 for instructions on preparing and running your RT experiment.

## Optimization strategies and troubleshooting

Refer to the pop-ups below for guidelines to optimize and troubleshoot your RT reaction.

RNA Sample Prep

RT Guidelines



Troubleshooting

Limited Warranty, Disclaimer,  
and Licensing Information

## Optional: Remove genomic DNA



The following example procedure describes the required to remove the genomic DNA (gDNA) from your samples for the RT and No RT Control reactions using the ezDNase™ Enzyme (supplied with Cat. No. 18091150 and 18091300; also available separately at [thermofisher.com](https://www.thermofisher.com)).

The volumes given are for a single 10-µL gDNA digestion reaction . For multiple reactions, prepare a master mix of components common to all reactions to minimize pipetting error, then dispense appropriate volumes into each reaction.

Step		Action										
1	 <p><b>Prepare gDNA digestion reaction mix (on ice)</b></p>	<p>For each RT reaction and No RT Control reaction, prepare a 10-µL gDNA reaction mix in an RNase-free tube on ice with the following components:</p> <table border="1"> <thead> <tr> <th>Component</th> <th>Volume</th> </tr> </thead> <tbody> <tr> <td>10× ezDNase™ Buffer</td> <td>1 µL</td> </tr> <tr> <td>ezDNase™ Enzyme</td> <td>1 µL</td> </tr> <tr> <td>Template RNA (1 pg to 2.5 µg total RNA)</td> <td>up to 8 µL</td> </tr> <tr> <td>Nuclease-free water</td> <td>to 10 µL</td> </tr> </tbody> </table>	Component	Volume	10× ezDNase™ Buffer	1 µL	ezDNase™ Enzyme	1 µL	Template RNA (1 pg to 2.5 µg total RNA)	up to 8 µL	Nuclease-free water	to 10 µL
Component	Volume											
10× ezDNase™ Buffer	1 µL											
ezDNase™ Enzyme	1 µL											
Template RNA (1 pg to 2.5 µg total RNA)	up to 8 µL											
Nuclease-free water	to 10 µL											
2	 <p><b>Digest gDNA</b></p>	Gently mix and incubate at 37°C for 2 minutes.										
3	 <p><b>Inactivate ezDNase™ enzyme</b></p>	If the RNA sample is to be used for RT-PCR of fragments ≥3 kb, incubate the sample for 5 minutes at 55°C in the presence of 10 mM DTT to inactivate the enzyme.										
4	 <p><b>Place the reaction on ice</b></p>	Chill the tube on ice to bring the sample to room temperature, then briefly centrifuge and place the tube on ice.										



## SuperScript™ IV First-Strand cDNA Synthesis Reaction

The following example procedure shows appropriate volumes for a single 20- $\mu\text{L}$  reverse transcription reaction. For multiple reactions, prepare a master mix of components common to all reactions to minimize pipetting error, then dispense appropriate volumes into each reaction tube prior to adding annealed template RNA and primers.

Step	Action										
<p><b>1</b></p>  <p><b>Anneal primer to template RNA</b></p>	<p>a. Combine the following components in a PCR reaction tube.  <b>Note:</b> Consider the volumes for all components listed in steps 1 and 2 to determine the correct amount of water required to reach your final reaction volume.</p> <table border="1" data-bbox="835 329 1858 548"> <thead> <tr> <th>Component</th> <th>Volume</th> </tr> </thead> <tbody> <tr> <td>50 <math>\mu\text{M}</math> Oligo d(T)<sub>20</sub> primer, 50 ng/<math>\mu\text{L}</math> random hexamers, or 2 <math>\mu\text{M}</math> gene-specific reverse primer</td> <td>1 <math>\mu\text{L}</math></td> </tr> <tr> <td>10 mM dNTP mix (10 mM each)</td> <td>1 <math>\mu\text{L}</math></td> </tr> <tr> <td>Template RNA (10 pg–5 <math>\mu\text{g}</math> total RNA or 10 pg–500 ng mRNA)</td> <td>up to 11 <math>\mu\text{L}</math></td> </tr> <tr> <td>DEPC-treated water</td> <td>to 13 <math>\mu\text{L}</math></td> </tr> </tbody> </table> <p>b. Mix and briefly centrifuge the components.            c. Heat the RNA-primer mix at 65°C for 5 minutes, and then incubate on ice for at least 1 minute.</p>	Component	Volume	50 $\mu\text{M}$ Oligo d(T) <sub>20</sub> primer, 50 ng/ $\mu\text{L}$ random hexamers, or 2 $\mu\text{M}$ gene-specific reverse primer	1 $\mu\text{L}$	10 mM dNTP mix (10 mM each)	1 $\mu\text{L}$	Template RNA (10 pg–5 $\mu\text{g}$ total RNA or 10 pg–500 ng mRNA)	up to 11 $\mu\text{L}$	DEPC-treated water	to 13 $\mu\text{L}$
Component	Volume										
50 $\mu\text{M}$ Oligo d(T) <sub>20</sub> primer, 50 ng/ $\mu\text{L}$ random hexamers, or 2 $\mu\text{M}$ gene-specific reverse primer	1 $\mu\text{L}$										
10 mM dNTP mix (10 mM each)	1 $\mu\text{L}$										
Template RNA (10 pg–5 $\mu\text{g}$ total RNA or 10 pg–500 ng mRNA)	up to 11 $\mu\text{L}$										
DEPC-treated water	to 13 $\mu\text{L}$										
<p><b>2</b></p>  <p><b>Prepare RT reaction mix</b></p>	<p>a. Vortex and briefly centrifuge the 5<math>\times</math> SSIV Buffer.            b. Combine the following components in a reaction tube.</p> <table border="1" data-bbox="835 724 1717 914"> <thead> <tr> <th>Component</th> <th>Volume</th> </tr> </thead> <tbody> <tr> <td>5<math>\times</math> SSIV Buffer</td> <td>4 <math>\mu\text{L}</math></td> </tr> <tr> <td>100 mM DTT</td> <td>1 <math>\mu\text{L}</math></td> </tr> <tr> <td>Ribonuclease Inhibitor</td> <td>1 <math>\mu\text{L}</math></td> </tr> <tr> <td>SuperScript™ IV Reverse Transcriptase (200 U/<math>\mu\text{L}</math>)</td> <td>1 <math>\mu\text{L}</math></td> </tr> </tbody> </table> <p>c. Cap the tube, mix, and then briefly centrifuge the contents.</p>	Component	Volume	5 $\times$ SSIV Buffer	4 $\mu\text{L}$	100 mM DTT	1 $\mu\text{L}$	Ribonuclease Inhibitor	1 $\mu\text{L}$	SuperScript™ IV Reverse Transcriptase (200 U/ $\mu\text{L}$ )	1 $\mu\text{L}$
Component	Volume										
5 $\times$ SSIV Buffer	4 $\mu\text{L}$										
100 mM DTT	1 $\mu\text{L}$										
Ribonuclease Inhibitor	1 $\mu\text{L}$										
SuperScript™ IV Reverse Transcriptase (200 U/ $\mu\text{L}$ )	1 $\mu\text{L}$										
<p><b>3</b></p>  <p><b>Combine annealed RNA and RT reaction mix</b></p>	<p>Add RT reaction mix to the annealed RNA.</p>										
<p><b>4</b></p>  <p><b>Incubate reactions</b></p>	<p>a. If using random hexamer, incubate the combined reaction mixture at 23°C for 10 minutes, and then proceed to step b.            If using oligo d(T)<sub>20</sub> or gene-specific primer, directly proceed to step b.            b. Incubate the combined reaction mixture at 50–55°C for 10 minutes.            c. Inactivate the reaction by incubating it at 80°C for 10 minutes.</p>										
<p><b>5</b></p>  <p><b>(Optional) Remove RNA</b></p>	<p><b>Note:</b> Amplification of some PCR targets (&gt;1 kb) may require removal of RNA. To remove RNA, add 1 <math>\mu\text{L}</math> <i>E. coli</i> RNase H, and incubate 37°C for 20 minutes.</p>										
<p><b>6</b></p>  <p><b>PCR amplification</b></p>	<p>Use your RT reaction immediately for PCR amplification or store it at –20°C.  <b>Note:</b> As a recommended starting point for PCR, reverse transcription reaction (cDNA) should compose 10% of the total reaction volume</p>										

## SuperScript™ IV Control Reactions - cDNA synthesis reaction


Follow the procedure below to perform the cDNA synthesis step of the SuperScript™ IV RT-PCR control reactions.

Step	Action	Action										
<p><b>1</b></p> 	<p><b>Anneal primer to template RNA</b></p>	<p>a. Prepare two tubes for annealing primer to template RNA. In each tube, combine the following:</p> <table border="1" data-bbox="835 240 1856 483"> <thead> <tr> <th>Component</th> <th>Volume</th> </tr> </thead> <tbody> <tr> <td>50 μM Oligo d(T)<sub>20</sub> primer</td> <td>1 μL</td> </tr> <tr> <td>10 mM dNTP mix (10 mM each)</td> <td>1 μL</td> </tr> <tr> <td>10 ng/μL total HeLa RNA (10 ng total)</td> <td>1 μL</td> </tr> <tr> <td>DEPC-treated water</td> <td>10 μL</td> </tr> </tbody> </table> <p>b. Mix and briefly centrifuge the components. c. Heat the RNA-primer mix at 65°C for 5 minutes, and then incubate on ice for at least 1 minute.</p>	Component	Volume	50 μM Oligo d(T) <sub>20</sub> primer	1 μL	10 mM dNTP mix (10 mM each)	1 μL	10 ng/μL total HeLa RNA (10 ng total)	1 μL	DEPC-treated water	10 μL
Component	Volume											
50 μM Oligo d(T) <sub>20</sub> primer	1 μL											
10 mM dNTP mix (10 mM each)	1 μL											
10 ng/μL total HeLa RNA (10 ng total)	1 μL											
DEPC-treated water	10 μL											
<p><b>2</b></p> 	<p><b>Prepare RT reaction mix</b></p>	<p>a. Vortex and briefly centrifuge the 5× SSIV Buffer. b. Prepare two reactions. In each reaction tube, combine the following:</p> <table border="1" data-bbox="835 675 1850 935"> <thead> <tr> <th>Component</th> <th>Volume</th> </tr> </thead> <tbody> <tr> <td>5× SSIV Buffer</td> <td>4 μL</td> </tr> <tr> <td>100 mM DTT</td> <td>1 μL</td> </tr> <tr> <td>Ribonuclease Inhibitor</td> <td>1 μL</td> </tr> <tr> <td>SuperScript™ IV Reverse Transcriptase (positive control) or DEPC-treated water (no RT control)</td> <td>1 μL</td> </tr> </tbody> </table> <p>c. Cap the tube, mix, and then briefly centrifuge the contents.</p>	Component	Volume	5× SSIV Buffer	4 μL	100 mM DTT	1 μL	Ribonuclease Inhibitor	1 μL	SuperScript™ IV Reverse Transcriptase (positive control) or DEPC-treated water (no RT control)	1 μL
Component	Volume											
5× SSIV Buffer	4 μL											
100 mM DTT	1 μL											
Ribonuclease Inhibitor	1 μL											
SuperScript™ IV Reverse Transcriptase (positive control) or DEPC-treated water (no RT control)	1 μL											
<p><b>3</b></p> 	<p><b>Combine annealed RNA and RT reaction mix</b></p>	<p>Add RT reaction mix to the annealed RNA.</p>										
<p><b>4</b></p> 	<p><b>Incubate reactions</b></p>	<p>a. Incubate the combined reaction mixture at 50°C for 10 minutes. b. Inactivate the reaction by incubating it at 80°C for 10 minutes.</p>										
<p><b>5</b></p> 	<p><b>Remove RNA</b></p>	<p>a. Add 1 μL <i>E. coli</i> RNase H and incubate 37°C for 20 minutes. b. Proceed to PCR amplification (page 5)</p>										

## SuperScript™ IV Control Reactions - PCR amplification

Follow the procedure below to perform the PCR amplification step of the SuperScript™ IV RT-PCR control reactions.

The example procedure below shows the appropriate volumes for a single 50- $\mu$ L reaction. For multiple reactions, prepare a master mix of components common to all reactions to minimize pipetting error, then dispense the appropriate volumes into each 0.2-mL or 0.5-mL PCR tube before adding template DNA and primers.

Step	Action																		
<p><b>1</b></p> 	<p><b>Assemble PCR amplification mix</b></p> <p>a. Prepare two reactions. In each tube, combine the following:</p> <table border="1" data-bbox="835 342 1850 857"> <thead> <tr> <th>Component</th> <th>Volume</th> </tr> </thead> <tbody> <tr> <td>DEPC-treated water</td> <td>37.8 <math>\mu</math>L</td> </tr> <tr> <td>10<math>\times</math> High Fidelity PCR Buffer</td> <td>5 <math>\mu</math>L</td> </tr> <tr> <td>50 mM MgSO<sub>4</sub></td> <td>2 <math>\mu</math>L</td> </tr> <tr> <td>10 mM dNTP mix (10 mM each)</td> <td>1 <math>\mu</math>L</td> </tr> <tr> <td>Control sense primer (10 <math>\mu</math>M) (5'-GCTCGTCGTCGACAACGGCTC-3')</td> <td>1 <math>\mu</math>L</td> </tr> <tr> <td>Control antisense primer (10 <math>\mu</math>M) (5'-CAAACATGATCTGGGTCATCTTCTC-3')</td> <td>1 <math>\mu</math>L</td> </tr> <tr> <td>cDNA from positive control reaction (step 5, page 4) or DEPC-treated water for no RT control</td> <td>2 <math>\mu</math>L</td> </tr> <tr> <td>Platinum™ Taq DNA Polymerase High Fidelity (5 U/<math>\mu</math>L)</td> <td>0.2 <math>\mu</math>L</td> </tr> </tbody> </table> <p>b. Mix gently by pipetting up and down and briefly centrifuge the components.</p>	Component	Volume	DEPC-treated water	37.8 $\mu$ L	10 $\times$ High Fidelity PCR Buffer	5 $\mu$ L	50 mM MgSO <sub>4</sub>	2 $\mu$ L	10 mM dNTP mix (10 mM each)	1 $\mu$ L	Control sense primer (10 $\mu$ M) (5'-GCTCGTCGTCGACAACGGCTC-3')	1 $\mu$ L	Control antisense primer (10 $\mu$ M) (5'-CAAACATGATCTGGGTCATCTTCTC-3')	1 $\mu$ L	cDNA from positive control reaction (step 5, page 4) or DEPC-treated water for no RT control	2 $\mu$ L	Platinum™ Taq DNA Polymerase High Fidelity (5 U/ $\mu$ L)	0.2 $\mu$ L
Component	Volume																		
DEPC-treated water	37.8 $\mu$ L																		
10 $\times$ High Fidelity PCR Buffer	5 $\mu$ L																		
50 mM MgSO <sub>4</sub>	2 $\mu$ L																		
10 mM dNTP mix (10 mM each)	1 $\mu$ L																		
Control sense primer (10 $\mu$ M) (5'-GCTCGTCGTCGACAACGGCTC-3')	1 $\mu$ L																		
Control antisense primer (10 $\mu$ M) (5'-CAAACATGATCTGGGTCATCTTCTC-3')	1 $\mu$ L																		
cDNA from positive control reaction (step 5, page 4) or DEPC-treated water for no RT control	2 $\mu$ L																		
Platinum™ Taq DNA Polymerase High Fidelity (5 U/ $\mu$ L)	0.2 $\mu$ L																		
<p><b>2</b></p> 	<p><b>Incubate reactions in a thermal cycler</b></p> <p>a. Place reaction mixture in preheated (94°C) thermal cycler.</p> <p>b. Perform PCR amplification using the following cycling parameters:</p> <table border="1" data-bbox="835 1032 1709 1318"> <thead> <tr> <th>Step</th> <th>Temperature</th> <th>Time</th> </tr> </thead> <tbody> <tr> <td>Initial denaturation</td> <td>94°C</td> <td>2 minutes</td> </tr> <tr> <td rowspan="3">35 PCR cycles</td> <td>Denature</td> <td>94°C</td> </tr> <tr> <td>Anneal</td> <td>55°C</td> </tr> <tr> <td>Extend</td> <td>68°C</td> </tr> <tr> <td>Hold</td> <td>4°C</td> <td>hold</td> </tr> </tbody> </table>	Step	Temperature	Time	Initial denaturation	94°C	2 minutes	35 PCR cycles	Denature	94°C	Anneal	55°C	Extend	68°C	Hold	4°C	hold		
Step	Temperature	Time																	
Initial denaturation	94°C	2 minutes																	
35 PCR cycles	Denature	94°C																	
	Anneal	55°C																	
	Extend	68°C																	
Hold	4°C	hold																	
<p><b>3</b></p> 	<p><b>Analyze with gel electrophoresis</b></p> <p>Analyze 10 <math>\mu</math>L of each reaction using agarose gel electrophoresis and ethidium bromide staining. A 353-bp band should be visible for the positive control reaction with RT. For the no RT control reaction, the same band should be <math>\leq</math> 50% in intensity when compared to the positive control.</p>																		