

BigDye™ Terminator v3.1 Matrix Standard Kit

SeqStudio™, 3500, and 3130 series instruments

Catalog Number 4336974

Pub. No. 4363117 Rev. C

WARNING! Read the Safety Data Sheets (SDSs) and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves. Safety Data Sheets (SDSs) are available from thermofisher.com/support.

Note: For safety and biohazard guidelines, see the “Safety” appendix in the the instrument user guide. Read the Safety Data Sheets (SDSs) and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves.

Product description

The BigDye™ Terminator v3.1 Matrix Standard Kit is used to perform spectral calibrations. The matrix standard contains four sizes of DNA fragments, each size labeled with a different fluorescent dye. The matrix standard can be used to perform spectral calibrations on the following instruments:

- SeqStudio™ Genetic Analyzer
- 3500/3500xL Genetic Analyzer
- 3130/3130xl Genetic Analyzer

Contents and storage

| Contents | Amount | Storage |
|---|--------|--|
| BigDye™ Terminator v3.1 Matrix Standard | 1 tube | Store at 2–8°C, protected from light. ^[1] Do not freeze. |

^[1] The kit is stable for six months when stored at 2–8°C.

Guidelines for use

- For more information on the use of matrix standards, see the instrument user guide or *DNA Fragment Analysis by Capillary Electrophoresis User Guide* (Pub. No. 4474504).
- **IMPORTANT!** Thoroughly mix the contents of the matrix standard tubes, then briefly centrifuge before use.
- To prepare the matrix standard dilution, combine the appropriate volumes of matrix standard and Hi-Di™ Formamide (Cat. No. 4311320). Dilution volumes vary depending on the specific application and instrument.
- Do not prepare the matrix standard more than 2 hours in advance.
- Do not add size standard to the matrix standard.
- **IMPORTANT!** Discard any unused reagent that has been diluted in Hi-Di™ Formamide.

Prepare the standard for the SeqStudio™ Genetic Analyzer

1. Combine the components for the array:

| Component | Volume |
|---------------------|-------------------|
| | 4-capillary array |
| Matrix standard | 1 µL |
| Hi-Di™ Formamide | 49 µL |
| Total volume | 50 µL |

2. Mix thoroughly, then centrifuge to bring the mixture to the bottom of the tube and eliminate air bubbles.
3. To denature the DNA fragments, incubate the mixture at 95°C for 2 minutes. Immediately place the mixture on ice.
4. Add 10 µL of the mixture to 4 wells (A1–D1) wells in a 96-well plate.
5. Cover the plate with a 96-well septa (Cat. No. 4315933), then centrifuge to bring the mixture to the bottom of the tube and eliminate air bubbles.
6. See the instrument user guide for specifics on setting up the run.

Prepare the standard for the 3500/3500xL Genetic Analyzer

1. Combine the components for the appropriate capillary array:

| Component | Volume | |
|---------------------|-------------------|--------------------|
| | 8-capillary array | 24-capillary array |
| Matrix standard | 2 µL | 5 µL |
| Hi-Di™ Formamide | 98 µL | 245 µL |
| Total volume | 100 µL | 250 µL |

2. Mix thoroughly, then centrifuge to bring the mixture to the bottom of the tube and eliminate air bubbles.
3. To denature the DNA fragments, incubate the mixture at 95°C for 2 minutes. Immediately place the mixture on ice.
4. Add 10 µL of the mixture to the appropriate number of wells in a 96-well plate:
 - **8-capillary array**—8 wells (A1–H1)
 - **24-capillary array**—24 wells (A1–H3)

Note: If you place the standard in other wells, specify the starting well in the software.
5. Cover the plate with a 96-well septa (Cat. No. 4315933), then centrifuge to bring the mixture to the bottom of the tube and eliminate air bubbles.
6. Assemble the plate with the Retainer and Base, then load on the instrument.
7. See the instrument user guide for specifics on setting up the run.

Prepare the standard for the 3130/3130x/ Genetic Analyzer

- Combine the components for the appropriate capillary array.

Table 1 16-capillary array

| Component | Volume | | |
|---------------------|---------------|---------------|---------------|
| | 36-cm array | 50-cm array | 80-cm array |
| Matrix standard | 10 µL | 5 µL | 10 µL |
| Hi-Di™ Formamide | 190 µL | 195 µL | 190 µL |
| Total volume | 200 µL | 200 µL | 200 µL |

Table 2 4-capillary array

| Component | Volume | | |
|---------------------|--------------|---------------|--------------|
| | 36-cm array | 50-cm array | 80-cm array |
| Matrix standard | 2.5 µL | 2.5 µL | 2.5 µL |
| Hi-Di™ Formamide | 47.5 µL | 97.5 µL | 47.5 µL |
| Total volume | 50 µL | 100 µL | 50 µL |

- Mix thoroughly, then centrifuge to bring the mixture to the bottom of the tube and eliminate air bubbles.
- To denature the DNA fragments, incubate the mixture at 95°C for 2 minutes. Immediately place the mixture on ice.
- Add 10 µL of the mixture to the appropriate number of wells in a 96-well plate.
 - 16-capillary array—16 wells (A1–H2)
 - 4-capillary array—4 wells (A1–D1)
- Assemble the plate with the Retainer and Base, then load on the instrument.
- See the instrument user guide for specifics on setting up the run.

Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale at www.thermofisher.com/us/en/home/global/terms-and-conditions.html. If you have any questions, please contact Life Technologies at www.thermofisher.com/support.



Thermo Fisher Scientific Baltics UAB | V.A. Graiciuno 8, LT-02241 | Vilnius, Lithuania

For descriptions of symbols on product labels or product documents, go to thermofisher.com/symbols-definition.

The information in this guide is subject to change without notice.

DISCLAIMER: TO THE EXTENT ALLOWED BY LAW, THERMO FISHER SCIENTIFIC INC. AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE, OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

Revision history: Pub. No. 4363117

| Revision | Date | Description |
|----------|------------------|---|
| C | 05 November 2018 | <ul style="list-style-type: none"> Update instrumentation Update licensing, trademarks, general style and format. |
| B | 18 August 2009 | Baseline for this revision |

Important Licensing Information: These products may be covered by one or more Limited Use Label Licenses. By use of these products, you accept the terms and conditions of all applicable Limited Use Label Licenses.

©2018 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.