

Safety data sheet
according to 1907/2006/EC, Article 31

Printing date 21.01.2021

Version number 4

Revision: 21.01.2021

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1 Product identifier

- **Trade name: 2-Mercaptoethanol**
- **Synonyma Thioglycol**



· **Article number:** 28625

· **CAS Number:**
60-24-2

· **EC number:**
200-464-6

· **Registration number** 01-2119517582-41

1.2 Relevant identified uses of the substance or mixture and uses advised against

- **Process category PROC15** Use as laboratory reagent
- **Application of the substance / the mixture** Laboratory chemicals

1.3 Details of the supplier of the safety data sheet

· **Manufacturer/Supplier:**
SERVA Electrophoresis GmbH
Carl-Benz-Str. 7
D-69115 Heidelberg
Tel.: +49/6221/13840-0
FAX: +49/6221/13840-10
msds.info@serva.de

· **Information department:** Product Safety department Tel.: +49 6221 13840-35

1.4 Emergency telephone number:

Medical Emergency Information in case of poisoning:
Poison Information Center Mainz - Phone: +49 (0) 6131 19240
(advisory service in German or English language)

SECTION 2: Hazards identification

2.1 Classification of the substance or mixture

· **Classification according to Regulation (EC) No 1272/2008**



GHS06

Acute Tox. 3 H301 Toxic if swallowed.
Acute Tox. 2 H310 Fatal in contact with skin.
Acute Tox. 3 H331 Toxic if inhaled.



GHS08

STOT RE 2 H373 May cause damage to organs through prolonged or repeated exposure.



GHS05

Eye Dam. 1 H318 Causes serious eye damage.



GHS09

Aquatic Acute 1 H400 Very toxic to aquatic life.
Aquatic Chronic 2 H411 Toxic to aquatic life with long lasting effects.

(Contd. on page 2)

Safety data sheet

according to 1907/2006/EC, Article 31

Printing date 21.01.2021

Version number 4

Revision: 21.01.2021

Trade name: 2-Mercaptoethanol

(Contd. of page 1)



GHS07

Skin Irrit. 2 H315 Causes skin irritation.
 Skin Sens. 1 H317 May cause an allergic skin reaction.

· **2.2 Label elements**

· **Labelling according to Regulation (EC) No 1272/2008**

The substance is classified and labelled according to the CLP regulation.

· **Hazard pictograms** GHS05, GHS06, GHS08, GHS09

· **Signal word** Danger

· **Hazard statements**

H301+H331 Toxic if swallowed or if inhaled.

H310 Fatal in contact with skin.

H315 Causes skin irritation.

H318 Causes serious eye damage.

H317 May cause an allergic skin reaction.

H373 May cause damage to organs through prolonged or repeated exposure.

H400 Very toxic to aquatic life.

H411 Toxic to aquatic life with long lasting effects.

· **Precautionary statements**

P273 Avoid release to the environment.

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P301+P310 IF SWALLOWED: Immediately call a POISON CENTER/ doctor.

P302+P352 IF ON SKIN: Wash with plenty of soap and water.

P304+P340 IF INHALED: Remove person to fresh air and keep comfortable for breathing.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

· **Labelling of packages where the contents do not exceed 125 ml**

· **Hazard pictograms** GHS05, GHS06, GHS08, GHS09

· **Signal word** Danger

· **Hazard statements**

H301+H331 Toxic if swallowed or if inhaled.

H310 Fatal in contact with skin.

H317 May cause an allergic skin reaction.

H373 May cause damage to organs through prolonged or repeated exposure.

· **Precautionary statements**

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P301+P310 IF SWALLOWED: Immediately call a POISON CENTER/ doctor.

P302+P352 IF ON SKIN: Wash with plenty of soap and water.

P304+P340 IF INHALED: Remove person to fresh air and keep comfortable for breathing.

· **2.3 Other hazards**

· **Results of PBT and vPvB assessment**

· **PBT:** PBT - assessment not available.

· **vPvB:** vPvB - assessment not available.

SECTION 3: Composition/information on ingredients

· **3.1 Chemical characterisation: Substances**

· **CAS No. Description:**

60-24-2 2-Mercaptoethanol

· **Identification number(s):**

· **EC number:** 200-464-6

· **Impurities and stabilising additives:**

· **Empirical formula:** C₂H₆O S

(Contd. on page 3)

GB

Safety data sheet
according to 1907/2006/EC, Article 31

Printing date 21.01.2021

Version number 4

Revision: 21.01.2021

Trade name: 2-Mercaptoethanol

(Contd. of page 2)

· **MW:** 78.13**SECTION 4: First aid measures**· **4.1 Description of first aid measures**· **General information**

Symptoms of poisoning may even occur after several hours; therefore medical observation for at least 48 hours after the accident.

Remove contaminated clothing.

Remove breathing apparatus only after contaminated clothing have been completely removed.

· **After inhalation** Supply fresh air or oxygen; call for doctor.

· **After skin contact**

Immediate wash with copious amounts of water and soap; rinse thoroughly; seek medical advice.

· **After eye contact**

Rinse opened eye for several minutes under running water. Remove present contact lenses, if easy to do, and continue rinsing. Consult ophthalmologist immediately.

· **After swallowing**

Wash out mouth. Drink plenty of water and supply fresh air.

Call for doctor immediately.

· **4.2 Most important symptoms and effects, both acute and delayed**

No further relevant information available.

· **4.3 Indication of any immediate medical attention and special treatment needed**

No further relevant information available.

SECTION 5: Firefighting measures· **5.1 Extinguishing media**· **Suitable extinguishing agents**

CO₂, extinguishing powder or water spray. Fight larger fires with water spray or alcohol resistant foam.

· **5.2 Special hazards arising from the substance or mixture**

In case of fire, the following can be formed, but not limited to:

Sulphur oxides (SO_x)

Hydrogen sulfide

Carbon monoxide and carbon dioxide

· **5.3 Advice for firefighters**

· **Protective equipment:** Wear self-contained respiratory protective device.

· **Additional information**

Collect contaminated fire fighting water separately. It must not enter the sewage system.

SECTION 6: Accidental release measures· **6.1 Personal precautions, protective equipment and emergency procedures**

Ensure adequate ventilation

Keep away from ignition sources

Avoid contact with the eyes and skin.

· **6.2 Environmental precautions:**

Inform respective authorities in case of seepage into water course or sewage system.

Do not allow to enter sewers/ surface or ground water.

· **6.3 Methods and material for containment and cleaning up:**

Dispose contaminated material as waste according to item 13.

Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).

· **6.4 Reference to other sections**

See Section 7 for information on safe handling

See Section 8 for information on personal protection equipment.

(Contd. on page 4)

Safety data sheet
according to 1907/2006/EC, Article 31

Printing date 21.01.2021

Version number 4

Revision: 21.01.2021

Trade name: 2-Mercaptoethanol

(Contd. of page 3)

See Section 13 for disposal information.

SECTION 7: Handling and storage

- **7.1 Precautions for safe handling**
Work only in fume cupboard.
Ensure good ventilation/exhaustion at the workplace.
Open and handle receptacle with care.
- **Information about protection against explosions and fires:** Keep respiratory protective device available.
- **7.2 Conditions for safe storage, including any incompatibilities**
- **Storage**
- **Requirements to be met by storerooms and receptacles:**
Store only in unopened original receptacles.
Store at +2 to +8 °C
- **Information about storage in one common storage facility:** Store away from oxidising agents.
- **Further information about storage conditions:**
Store under lock and key and with access restricted to technical experts or their assistants only.
Keep receptacle tightly sealed and store in dry conditions.
- **7.3 Specific end use(s)** No further relevant information available.

SECTION 8: Exposure controls/personal protection

- **8.1 Control parameters**
- **Additional information about design of technical systems:** No further data; see item 7.
- **Components with limit values that require monitoring at the workplace:** Not required.
- **DNELs**
worker: long-term-systemic effects, dermal: 0,6 mg/kg
worker: long-term-systemic effects, inhalation: 4 mg/m³
- **PNECs**
PNEC fresh water: 0,0004 mg/l
PNEC fresh water sediment: 0,00084 mg/kg
PNEC marine water: 0,00004 mg/l
PNEC Soil: 0,29175 mg/kg
- **Additional information:** The lists that were valid during the creation were used as basis.
- **8.2 Exposure controls**
- **Personal protective equipment**
- **General protective and hygienic measures**
Keep away from foodstuffs, beverages and feed.
Immediately remove all soiled and contaminated clothing
Store protective clothing separately.
Avoid contact with the eyes and skin.
Wash hands before breaks and at the end of work.
Avoid contact with the skin.
- **Breathing equipment:**
Short term filter device:
Filter A/P3
- **Protection of hands:**
The glove material has to be impermeable and resistant to the product/ the substance/ the preparation.
Protective gloves.
Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/ the chemical mixture.
Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation

(Contd. on page 5)

Safety data sheet
according to 1907/2006/EC, Article 31

Printing date 21.01.2021

Version number 4

Revision: 21.01.2021

Trade name: 2-Mercaptoethanol

(Contd. of page 4)

- **Material of gloves**

The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer.

- **Penetration time of glove material**

The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.

- **For the permanent contact of a maximum of 15 minutes gloves made of the following materials are suitable:**

Nitrile rubber, NBR

PVC gloves

Butyl rubber, BR

- **Eye protection:** Tightly sealed goggles.

- **Body protection:** Protective work clothing.

SECTION 9: Physical and chemical properties

- **9.1 Information on basic physical and chemical properties**

- **General Information**

- **Appearance:**

Form:	Liquid
Colour:	Colourless
Odour:	Unpleasant
Odour threshold:	Not determined.

· **pH-value:** 5

- **Change in condition**

Melting point/freezing point: undetermined
Initial boiling point and boiling range: 155-160 °C

· **Flash point:** 74 °C

· **Flammability (solid, gaseous)** Not applicable.

· **Ignition temperature:** 295 °C

· **Decomposition temperature:** Not determined.

· **Self igniting:** Not determined.

· **Explosive properties:** Product does not present an explosion hazard.

- **Explosion limits:**

Lower: 2.3 Vol %
Upper: 18 Vol %

· **Vapour pressure at 20 °C:** 1.3 hPa

· **Density at 20 °C:** 1.12 g/cm³

· **Relative density** Not determined.

· **Vapour density** Not determined.

· **Evaporation rate** Not determined.

- **Solubility in / Miscibility with**

Water: Fully miscible

· **Partition coefficient: n-octanol/water:** Not determined.

- **Viscosity:**

dynamic at 20 °C: 3.22 mPas

kinematic: Not determined.

(Contd. on page 6)

Safety data sheet
according to 1907/2006/EC, Article 31

Printing date 21.01.2021

Version number 4

Revision: 21.01.2021

Trade name: 2-Mercaptoethanol

(Contd. of page 5)

· **9.2 Other information** *No further relevant information available.*

SECTION 10: Stability and reactivity

- **10.1 Reactivity** *No further relevant informations available*
- **10.2 Chemical stability**
- **Thermal decomposition / conditions to be avoided:**
No decomposition if used and stored according to specifications.
- **10.3 Possibility of hazardous reactions** *Vapours can form flammable and explosive mixtures with air.*
- **10.4 Conditions to avoid**
*Avoid high temperatures, flames, sparks
moisture*
- **10.5 Incompatible materials:** *Avoid contact with strong oxidizers and reducing agents.*
- **10.6 Hazardous decomposition products:** *In case of fire: See Section 5*

SECTION 11: Toxicological information

- **11.1 Information on toxicological effects**
 - **Acute toxicity**
*Toxic if swallowed or if inhaled.
Fatal in contact with skin.*
- | · LD/LC50 values that are relevant for classification: | | |
|---|----------|------------------------|
| Oral | LD50 | 133 mg/kg (rat) |
| Dermal | LD50 | 112-224 mg/kg (rabbit) |
| Inhalative | LC50/4h | 2.1 mg/l (rat) |
| | LC50/96h | 37 mg/l (Forelle) |
- **Primary irritant effect:**
 - **Skin corrosion/irritation**
Causes skin irritation.
 - **Serious eye damage/irritation**
Causes serious eye damage.
 - **Respiratory or skin sensitisation**
May cause an allergic skin reaction.
 - **Additional toxicological information:**
 - **CMR effects (carcinogenicity, mutagenicity and toxicity for reproduction)**
 - **Germ cell mutagenicity** *Based on available data, the classification criteria are not met.*
 - **Carcinogenicity** *Based on available data, the classification criteria are not met.*
 - **Reproductive toxicity** *Based on available data, the classification criteria are not met.*
 - **STOT-single exposure** *Based on available data, the classification criteria are not met.*
 - **STOT-repeated exposure**
May cause damage to organs through prolonged or repeated exposure.
 - **Aspiration hazard** *Based on available data, the classification criteria are not met.*

SECTION 12: Ecological information

- **12.1 Toxicity**
 - **Aquatic toxicity:**
- | | |
|----------|--|
| EC50/48h | 0.4 mg/l (<i>Daphnia magna</i>) |
| EC50/72h | 19 mg/l (<i>Scenedesmus subspicatus</i>) |
- **12.2 Persistence and degradability** *No further relevant information available.*
 - **Other information:** *Biodegradability: >70% in 28d (OECD 309)*

(Contd. on page 7)

Safety data sheet
according to 1907/2006/EC, Article 31

Printing date 21.01.2021

Version number 4

Revision: 21.01.2021

Trade name: 2-Mercaptoethanol



(Contd. of page 6)

- **12.3 Bioaccumulative potential**
bioaccumulation potential is not to be expected
log Pow = -0,056
- **12.4 Mobility in soil** No further relevant information available.
- **Ecotoxicological effects:**
- **Remark:** Very toxic to aquatic life with long lasting effects.
- **Additional ecological information:**
- **General notes:**
Do not allow product to reach ground water, water course or sewage system.
Water danger class 3 (German Regulation) (Assessment by list): extremely hazardous for water.
- **12.5 Results of PBT and vPvB assessment**
- **PBT:** PBT - assessment not available.
- **vPvB:** vPvB - assessment not available.
- **12.6 Other adverse effects** Toxicity to bacteria (EC50, 17 h): 125 mg/l (*Pseudomonas putida*)

SECTION 13: Disposal considerations

- **13.1 Waste treatment methods**
- **Recommendation**
Must be specially treated adhering to official regulations.
Must not be disposed of together with household garbage. Do not allow product to reach sewage system.
- **Uncleaned packagings:**
- **Recommendation:**
Disposal of uncleaned packagings must be made according to official regulations in the same manner as the product.
- **Recommended cleansing agent:** Water, if necessary with cleansing agents.

SECTION 14: Transport information

- | | |
|---|--|
| 14.1 UN-Number | UN2966 |
| ADR, IMDG, IATA | |
| 14.2 UN proper shipping name | 2966 THIOGLYCOL, ENVIRONMENTALLY HAZARDOUS |
| ADR | THIOGLYCOL, MARINE POLLUTANT |
| IMDG | THIOGLYCOL |
| IATA | |
| 14.3 Transport hazard class(es) | |
| ADR, IMDG | |
|  | |
| Class | 6.1 Toxic substances. |
| Label | 6.1 |
| IATA | |
|  | |
| Class | 6.1 Toxic substances. |
| Label | 6.1 |

(Contd. on page 8)

Safety data sheet
according to 1907/2006/EC, Article 31

Printing date 21.01.2021

Version number 4

Revision: 21.01.2021

Trade name: 2-Mercaptoethanol

(Contd. of page 7)

· 14.4 Packing group · ADR, IMDG, IATA	II
· 14.5 Environmental hazards: · Marine pollutant: · Special marking (ADR):	<i>Environmentally hazardous substance, liquid; Marine Pollutant</i> <i>Ja</i> <i>Symbol (fish and tree)</i> <i>Symbol (fish and tree)</i>
· 14.6 Special precautions for user · Hazard identification number (Kemler code): · EMS Number: · Stowage Category	<i>Warning: Toxic substances.</i> <i>60</i> <i>F-A,S-A</i> <i>A</i>
· 14.7 Transport in bulk according to Annex II of Marpol and the IBC Code	<i>Not applicable.</i>
· Transport/Additional information:	
· ADR · Limited quantities (LQ) · Excepted quantities (EQ)	<i>100 ml</i> <i>Code: E4</i> <i>Maximum net quantity per inner packaging: 1 ml</i> <i>Maximum net quantity per outer packaging: 500 ml</i>
· Transport category · Tunnel restriction code	<i>2</i> <i>D/E</i>
· IMDG · Limited quantities (LQ) · Excepted quantities (EQ)	<i>100 ml</i> <i>Code: E4</i> <i>Maximum net quantity per inner packaging: 1 ml</i> <i>Maximum net quantity per outer packaging: 500 ml</i>
· UN "Model Regulation":	<i>UN 2966 THIOGLYCOL, 6.1, II, ENVIRONMENTALLY HAZARDOUS</i>

SECTION 15: Regulatory information

- **15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture**
 - **Directive 2012/18/EU**
 - **Named dangerous substances - ANNEX I** Substance is not listed.
 - **Seveso category**
H2 ACUTE TOXIC
E1 Hazardous to the Aquatic Environment
 - **Qualifying quantity (tonnes) for the application of lower-tier requirements** 50 t
 - **Qualifying quantity (tonnes) for the application of upper-tier requirements** 200 t
 - **REGULATION (EC) No 1907/2006 ANNEX XVII** Conditions of restriction: 3
 - **DIRECTIVE 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic equipment – Annex II**
Substance is not listed.
 - **National regulations**
 - **Technical instructions (air):**
- | Class | Share in % |
|-------|------------|
| I | 80-100 |
- **Water hazard class:** Water danger class 3 (Assessment by list): extremely hazardous for water.

(Contd. on page 9)

GB

Safety data sheet
according to 1907/2006/EC, Article 31

Printing date 21.01.2021

Version number 4

Revision: 21.01.2021

Trade name: 2-Mercaptoethanol

(Contd. of page 8)

· **15.2 Chemical safety assessment:** A Chemical Safety Assessment has not been carried out.

SECTION 16: Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

· **Department issuing SDS:** Product safety department

· **Contact:** +49 6221 13840-35

· **Abbreviations and acronyms:**

RID: Règlement international concernant le transport des marchandises dangereuses par chemin de fer (Regulations Concerning the International Transport of Dangerous Goods by Rail)

ICAO: International Civil Aviation Organisation

PBT: persistent, bioaccumulative, toxic substance (REACH)

vPvB: very persistent, very bioaccumulative substance (REACH)

REACH: Regulation concerning the Registration, Evaluation, Authorisation and Restriction of Chemicals

CLP: Regulation on classification, labelling and packaging of substances and mixtures

bw: body weight

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

IATA: International Air Transport Association

GHS: Globally Harmonised System of Classification and Labelling of Chemicals

EINECS: European Inventory of Existing Commercial Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

DNEL: Derived No-Effect Level (REACH)

PNEC: Predicted No-Effect Concentration (REACH)

LC50: Lethal concentration, 50 percent

LD50: Lethal dose, 50 percent

PBT: Persistent, Bioaccumulative and Toxic

vPvB: very Persistent and very Bioaccumulative

Acute Tox. 3: Acute toxicity - oral – Category 3

Acute Tox. 2: Acute toxicity - dermal – Category 2

Skin Irrit. 2: Skin corrosion/irritation – Category 2

Eye Dam. 1: Serious eye damage/eye irritation – Category 1

Skin Sens. 1: Skin sensitisation – Category 1

STOT RE 2: Specific target organ toxicity (repeated exposure) – Category 2

Aquatic Acute 1: Hazardous to the aquatic environment - acute aquatic hazard – Category 1

Aquatic Chronic 2: Hazardous to the aquatic environment - long-term aquatic hazard – Category 2

GB