

## BD3 Cycle Sequencing Mix User Manual

### Product catalog numbers:

Catalog number	Product Name	Unit Size
6304010	BD3 Cycle Sequencing Mix	100µl
6305010	BD3 Cycle Sequencing Mix	125rxn
6305050	BD3 Cycle Sequencing Mix	125rxn*5
6305080	BD3 Cycle Sequencing Mix	125rxn*8
6305200	BD3 Cycle Sequencing Mix	125rxn*20
6306040	BD3 Cycle Sequencing Mix	2500rxn*2

### Product description:

BD3 Cycle Sequencing Mix is used to set up cycle sequencing reactions and contains all the components except for sequencing primers and DNA templates. The Mix is a high-performance alternative to BigDye® Terminator v3.1 Cycle Sequencing Mix.

### Product usage, performance, and storage conditions:

BD3 Cycle Sequencing Mix is used as a "drop in" product for your standard SOP for cycle sequencing reaction setup. There is no need to make changes in SOPs to achieve high-quality sequencing results using either downstream Onelife products or other products you normally use.

The product is shipped with dry ice and should be stored at -20°C. The product has a shelf life of 18 months when stored as suggested.