



NOT TO SCALE
120.00 ±0.25

99.00 ±0.25

Copyright Thermo Fisher Scientific
This drawing and/or any technical information relating to it is issued in strict confidence and may not be reproduced or used for manufacture without written permission. Upon receipt of which it is the users responsibility to seek confirmation of any dimensions and production tolerances.



VIAL FORM 96pos



TITLE
Thermo-Fast® 96,
Semi-Skirted PCR Plate

All dimensions in mm

DRG. NO.	ISSUE	DATE
AB-0900	F	7/2/13